


VARIO report

DICOM Viewer



The Intuitive Viewing Solution

- Access to all available images within seconds
- Easy, intuitive operation
- Seamless workflow integration with your existing practice management system



VARIO report - at a glance:

Easy, intuitive operation

Due to its flexible GUI, a modality-specific set of tools is provided to make the operation of the system easy and intuitive. The GUI can also be adapted to specific user requirements.

Modality-specific setting

For maximum flexibility, display modes can be configured to the requirements of different modalities.

Organ-specific setting

Organ-specific display configurations are possible to optimize the workflow in different medical scenarios.

Multiple Display

The simultaneous display of several examinations is supported.

Display of image information

Image information displays can be configured freely to show a selection of DICOM header tags.

Lookup Tables

To optimize image display, a choice of lookup options (modality LUT, VOI LUT, user-specific LUT) is available for image display and printing. All lookup options can be configured at runtime.

DICOM Print

In combination with DICOM Print, a substantial set of functions for printing is provided, including but not limited to annotations, print formats and media selection.

Paper Print

All DICOM Print options are available for paper printing as well. Special lookup tables can be implemented to compensate for gray scale interpretation of paper printers. All MS Windows-compatible paper printers are supported.

Patient-CD

Patient-CDs can be created quickly and easily. The SST VARIO report user interface allows for CD creation with only two mouse clicks. A JAVA-based DICOM viewer is part of the data set stored on the CD. This allows for displaying DICOM images without additional software.



Schaef Systemtechnik GmbH
Falkenstrasse 22
91580 Petersaurach
Germany

Phone.: +49 98 72 / 20 02
Fax: +49 98 72 / 28 23

eMail: info@sstmed.com
web: www.sstmed.com

Ask us about our other products:

- **VARIO store** - the high-speed DICOM PACS solution
- **ACC** - the professional image acquisition software
- **VARIO live** - the grabber software for ultrasound, endoscopy and arthroscopy